

## Hip Belt Extension: Universal

The ISO Preferred Lumbar Belt Extension: Universal provides an extension to the hip belt for waist sizes over 56 inches.



ISO Preferred Hip Belt Extension: Universal (HI204) Applies to two different types of Hip Braces, the ISO Preferred Post OP Hip Brace & The Clam Shell models.

### CLAM SHELL MODELS: - Belt attaches to brace with clamshell gripping material:

ISO PREFERRED  
DR MEDICAL

ISO-HI206  
HB-166D

**BEFORE APPLICATION:** You will be applying 2 Belt Extensions, one on each side of the back brace. Detach and extend one side of the hip brace belt at a time.

**WHEN COMPLETE, MAKE SURE THE FRONT SIDE WITH THE HOOK AND LOOP FASTENING IS ON THE RIGHT SIDE OF THE PATIENT'S BODY WHEN WORN.**

1. **APPLY THE EXTENSION:** On one side, open up the clam shell fastener, detach the hip belt and apply Extension Belt, facing down.
2. Attach the original hip belt to the Belt Extension by positioning leftward or rightward to add the desired number of inches for this side. Secure original belt to the Extension. Apply Belt Extension to the other side of the brace making sure both sides are even. Apply brace as instructed on product instructions.

1



2



# Hip Brace with Flexion/Extension/Abduction: Universal

The ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal: reduces the time it takes to heal discomfort from injuries or post-surgery. Please follow your orthopedic professional's instructions explicitly and adjust or discontinue use only as directed to achieve maximum benefit for your concerns.



ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal (HI206) fits right or left side and is intended to provide support for Hip Orthosis, Abduction/adduction control of hip joint, Post-Op hip abduction/adduction, Post-Op and Pre-Op support, Hip Dysplasia, Arthroplasty, Hip revisions and surgeries, Additional support for lower extremity, Pelvic control, Hip joint support with adjustable flexion/extension and abduction/adduction control.

1. The belt at waist is designed to fit sizes Small-2XL (optional extension available: ISO-HI204). Adjust waist to size by removing belt from each side of back panel at tab.



Select size, using the size markings on the belt. Fold over at correct size and reattach to back panel. (Cut as needed once proper fit is completed at Step 7). Attach leg support to right or left side as prescribed by your physician. To apply wrap belt around waist and using thumb pulls, secure to waist. Secure thumb pulls to belt.

2. Align hip attachment with patient's hip. Secure the hip attachment to waist belt using the hook and loop tabs.
3. For better compression and fit, pull string through the opening at the top of each plastic hip piece.
4. Adjust and position thigh supports by loosening, adjusting and tightening the thigh bolt as needed with the wrench provided.
5. To fit to patient hip and align leg at midline, loosen the bolt on the back side of ball joint and adjust abduction/adduction as needed.
6. Adjust Range of Motion (ROM) by loosening the small bolt at the top of ball joint assembly. After setting the extension and flexion tabs to patient's needs, tighten small bolt at the top to maintain settings.
7. Secure to patient legs with hook and loop straps. Straps may be shortened and cut as needed. Cut excess length from waist belt if desired, always leaving at least 3" overlap (or more) on each side to accommodate changes in waist size.



Adjust and shorten straps as needed to provide optimum support and comfort. **DO NOT OVERTIGHTEN.**

## ISO Hip Belt Extension: Universal (Optional-HI204)

For additional length at waist securely fasten the (optional) ISO Hip Belt Extension: Universal to each side of the ISO Hip Brace Belt.





## INTRODUCTION

Your **ISO Preferred Hip Belt Extension: Universal (HI204)** is designed to provide your body with the support and stability that it needs to successfully manage or eliminate your mobility or pain concerns. It is critical that you follow your orthopedic professional's directions closely to ensure that the apparatus is worn properly and in the manner that will be most beneficial in addressing your issues. For correct application, please read the enclosed instructions and watch the on-line instructional video (if available) at [www.ISOpreferred.com](http://www.ISOpreferred.com) so that this device can provide you with superior long-term usability and care.

## USAGE RECOMMENDATIONS

It is vital that you read the included instructions and when applicable watch the instructional video for the **ISO Preferred Hip Belt Extension: Universal (HI204)**. The brace is intended for application by qualified professionals only. It is not recommended for patients to fit or apply or adjust their new brace without an orthopedic professional's expertise, unless otherwise instructed. Further questions regarding adjustments or changes, frequency of use, and/or duration of use should be presented to your orthopedic professional, as they are familiar with the specific needs and conditions, for which this device was prescribed.

## CARE AND CLEANING INSTRUCTIONS

Please do not place your brace in the washing machine or dryer at any time. When preparing to clean your **ISO Preferred Hip Belt Extension: Universal (HI204)** you will want to separate the parts where applicable and care for them individually.

- Do not dry clean, iron, bleach, or use an abrasive cleaning instrument on any part of the brace. Failure to follow these instructions may damage your brace and void the manufacturer's warranty and it may no longer work properly.
- Hand wash each section in cold water with a mild detergent. Lay flat to air dry. Allow all parts to air dry prior to reassembling the device.

## CAUTION

The **ISO Preferred Hip Belt Extension: Universal (HI204)** is for single patient use only. While every effort has been made to obtain maximum strength, function and durability, there is no guarantee that injury will be prevented through the use of this device, especially if the instructions are not followed correctly.

**DO NOT OVERTIGHTEN.**



## INTRODUCTION

Your **ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal (HI206)** is designed to provide your body with the support and stability that it needs to successfully manage or eliminate your mobility concerns. It is critical that you follow your orthopedic professional's directions closely to ensure that the apparatus is worn properly and in the manner that will be most beneficial in addressing your issues. For correct application, please read the enclosed instructions and watch the on-line instructional video (if available) at [www.ISOPreferred.com](http://www.ISOPreferred.com) so that this device can provide you with superior long-term usability and care.

## USAGE RECOMMENDATIONS

It is vital that you read the included instructions and when applicable watch the instructional video for the **ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal (HI206)**. The brace is intended for application by qualified professionals only. It is not recommended for patients to fit or apply or adjust their new brace without an orthopedic professional's expertise, unless otherwise instructed. Further questions regarding adjustments or changes, frequency of use, and/or duration of use should be presented to your orthopedic professional, as they are familiar with the specific needs and conditions, for which this device was prescribed.

## CARE AND CLEANING INSTRUCTIONS

Please do not place your brace in the washing machine or dryer at any time. When preparing to clean your **ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal (HI206)** you will want to separate the parts where applicable and care for them individually.

- Do not dry clean, iron, bleach, or use an abrasive cleaning instrument on any part of the brace. Failure to follow these instructions may damage your brace and void the manufacturer's warranty and it may no longer work properly.
- Remove the inner padding
- Hand wash each section in cold water with a mild detergent. Lay flat to air dry.
- Allow all parts to air dry prior to reassembling the device.

## CAUTION

The **ISO Preferred Hip Brace with Flexion/Extension/Abduction: Universal (HI206)** is for single patient use only. While every effort has been made to obtain maximum strength, function and durability, there is no guarantee that injury will be prevented through the use of this device, especially if the instructions are not followed correctly. **DO NOT OVERTIGHTEN.**